NYLON

Instructions for Use

NYLON

(NONABSORBABLE SURGICAL SUTURES, USP / EP)

DESCRIPTION

Nylon is a nonabsorbable sterile surgical suture composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6,6. Nylon is monofilament sutures, clear or dyed black or blue. Nylon meets all requirements established by the United States Pharmacopeia (USP) and European Pharmacopeia (EP) for Nonabsorbable Surgical Suture.

INDICATIONS

Nylon is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

ACTIONS

Nylon elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While is not absorbed, progressive hydrolysis of the nylon in vivo may result in gradual loss of up to 20% of its tensile strength per year.

CONTRAINDICATIONS

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, Nylon should not be used where permanent retention of tensile strength is required, as in fixation of intraocular lenses or synthetic vascular grafts.

WARNINGS

Prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Do not resterilize. Discard open, unused sutures. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Nylon for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, gradual loss of tensile strength over time, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema and erythema at the wound site.

HOW SUPPLIED

Nylon is available in sizes 11-0 through 2. Nylon is undyed or dyed black with Logwood extract or blue with FD & C blue #2. The suture is supplied sterile in pre-cut lengths or ligating reels; both non-needled and affixed to various needle types; in one-dozen boxes.

Caution: National or local regulations may restrict this device to sale, distribution and use by or on the order of a physician, veterinarian or other licensed practitioner.

ProAdvantage by NDC
Made in BUSAN, KOREA

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