SILK

Instructions for Use

SILK

(NONABSORBABLE SURGICAL SUTURES (USP / EP))

DESCRIPTION

Silk is a nonabsorbable, sterile, surgical suture composed of an organic protein called fibroin. This protein is derived from the domesticated species *Bombyx mori* (B. mori) of the family Bombycidae.

Silk sutures are braided or twisted and available white or dyed (black or blue).

Silk meets all requirements established by the United States Pharmacopeia (USP) and European Pharmacopeia (EP) for Nonabsorbable Surgical Suture.

INDICATIONS

Silk is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

ACTIONS

Silk elicits an acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While Silk is not absorbed, progressive degradation of the proteinaceous silk fiber in vivo may result in gradual loss of all of the suture's tensile strength within one year.

CONTRAINDICATIONS

The use of this suture is contraindicated in patients with known sensitivities or allergies to silk.

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, Silk should not be used where permanent retention of tensile strength is required, as in fixation of vascular prosthesis.

WARNING

Prolonged contact of this or any other suture with salt solutions, such as those found in urinary or biliary tracts, may result in calculus formation. Do not resterilize. Discard open, unused sutures.

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Silk for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, gradual loss of all tensile strength over time, allergic response in patients that are known to be sensitive to silk, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, acute inflammatory tissue reaction, and pain, edema, and erythema at the wound site.

HOW SUPPLIED

Silk is available in sizes 10-0 through 2. Silk is white or dyed black with Logwood extract or blue with D & C Blue #9. The suture is supplied sterile in pre-cut lengths or ligating reels; both non-needled and affixed to various needle types; in one dozen boxes.

Caution: National or local regulations may restrict this device to sale, distribution and use by or on the order of a physician, veterinarian or other licensed practitioner.

ProAdvantage by NDC
Made in BUSAN, KOREA

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