



Blue, Nonwoven

CSR Wrap

(Central Supply Room Wrap)

A blue, protective wrap made from natural cellulose fibers for sterilizing products. Ideally suited for operating rooms, central supply, kit packing and other areas where equipment is wrapped prior to sterilization. Material is strong and tear resistant, yet soft and flexible for easy wrapping. Can be used with various methods of sterilization including Steam, ETO, gamma irradiation, and E-Beam. Provides the highest barrier to air and waterborne contamination.

- Autoclavable
- Disposable
- Latex-Free



Item #	Size	Qty
P011515	15" x 15"	100 sheets/bag
P012020	20" x 20"	100 sheets/bag
P012424	24" x 24"	100 sheets/bag

Distributed By: