

# Glutaraldehyde

## Sterilizing & Disinfecting Solution

### 14-Day

- Ideal for lensed instruments such as flexible and/or rigid scopes
- High level disinfection is achieved in 45 minutes
- Sterilization is achieved in 10 hours
- Compatible with heat-sensitive medical equipment and instrumentation
- Contains no surfactants
- No dilution required

Item #	Description	UOM
<b>N099001</b>	<b>14-Day Solution, 2.6% Glutaraldehyde, High Level Disinfectant/Sterilant, Gallon</b>	<b>4/cs</b>

### 28-Day

- Recommended for immersible medical and heat sensitive equipment and instruments
- High level disinfection is achieved in 90 minutes
- Sterilization is achieved in 10 hours
- Contains a surfactant for better cleaning
- Contains a rust-inhibitor to reduce corrosion
- No dilution required

Item #	Description	UOM
<b>N099002</b>	<b>28-Day Solution, 2.5% Glutaraldehyde, High Level Disinfectant/Sterilant, Gallon</b>	<b>4/cs</b>



*For quality control, we recommend the use of MetriTest® Test Strips. Use MetriTest® 1.5% for use with 14-Day (#10-303). Use MetriTest® 1.8% for use with 28-Day (#10-304). Strips should be used before immersing instruments.*



**WARNING:** This product can expose you to chemicals including methanol, which is known to the State of California to cause birth defects or other reproductive harm.

**Distributed by:**