

Disposable

Ear Curettes

Pro Advantage double ended ear curettes allow healthcare professionals to utilize a single device for safe and effective removal of impacted or occluded cerumen from the ear canal.

- Unique double ended design features choice of paddle or loop end with different sized tips on each end
- “No Roll” triangular handles are easy to grip and allow for firm control of the tips
- Constructed of plastic for minimal trauma to the ear
- Orientation marks allow for alignment of tips from outside the ear
- Elongated tapers permit better view into the ear canal
- Handy dispenser box easily fits into drawer
- Single use only



Item #	Description	Quantity
P040002	Loop Ear Curette	50/box
P040004	Paddle Ear Curette	50/box