

Throat Swab Strep A Test Strips

**A^{Pro}
ADVANTAGE[®]**
by **NDC**

A rapid test for the qualitative detection of Strep A antigen in throat swab specimens.

Lateral flow immunoassay technology identifies antibodies specific for whole cell Lancefield Group A Streptococcus to selectively detect Strep A carbohydrate antigen in a throat swab specimen. Aids in the early diagnosis of Group A Streptococcal infection.

- 97% Sensitivity; 95% Specificity
- Accurate results in 5 minutes
- Quickly diagnose and treat infections to help avoid serious complications
- CLIA waived and FDA 510(k) cleared to market
- Room temperature storage



CPT Code: 87880

Item #	Description
P080033	Strep A Test-Throat Swab Kit <i>Includes: 25 test strips, 25 swabs, 25 extraction tubes, 1-10mL reagent A, 1-10mL reagent B, 1-1mL positive control, 1-1mL negative control, package insert, procedure card. 25 tests per box</i>

Distributed By: